

E 497

Nozzle

for optimum loading of double and single tube breathing masks. Height 195 mm.

EAN: 4002514811253 / Material number: 07759470 / Old Material Number: 69749701D



- H 195 mm, with pipe 109 mm at the side

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing anaesthetic instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	195
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	3 15/16 (100)
External dimensions, gross depth in in (mm)	5 15/16 (150)
Net weight in kg	1 (0.06)
Gross weight in kg	1 (0.07)