

E 900-5/2
Mobile unit
for the optimum cleaning of instruments on 5 levels.



- Loading capacity: 15 DIN mesh trays*
- The loading areas for the individual levels are different
- Magnetic strip for wash cart recognition (incl. magnets ML)

* DIN mesh tray 485 x 254 x 50 mm

EAN: 4002514816418 / Material number: 07765760 / Old Material Number: 69590006D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	15
SPRI mesh trays [number]	10
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	23 7/8 (605)
External dimensions, net width in in. (mm)	25 1/4 (640)
External dimensions, net depth in in. (mm)	31 1/8 (790)
External dimensions, gross height in in. (mm)	26 13/16 (680)
External dimensions, gross width in in. (mm)	22 1/16 (560)
External dimensions, gross depth in in. (mm)	25 5/8 (650)
Net weight in lb (kg)	29 (12.85)
Gross weight in lb (kg)	29 (12.85)
Accessory pack	
Magnetic strip for automatic program recognition (incl. magnets) [number]	1