

E 935/2

Mobile unit

for the optimum cleaning of instruments on 3 levels.



- Loading capacity: 9 DIN mesh trays*
- Equipped with 2 integrated spray arms
- The loading areas for the individual levels are different
- Magnetic strip for wash cart recognition (incl. magnets)

* DIN mesh tray 485 x 254 x 50 mm

EAN: 4002514816432 / Material number: 07765780 / Old Material Number: 69593503D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	9
SPRI mesh trays [number]	6
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	20 11/16 (524)
External dimensions, net width in in. (mm)	25 1/4 (640)
External dimensions, net depth in in. (mm)	31 1/8 (790)
External dimensions, gross height in in. (mm)	28 3/8 (720)
External dimensions, gross width in in. (mm)	32 11/16 (830)
External dimensions, gross depth in in. (mm)	30 3/4 (780)
Net weight in lb (kg)	59 (26.50)
Gross weight in lb (kg)	67 (30.00)
Accessory pack	
Magnetic strip for automatic program recognition (incl. magnets) [number]	1