

E 975/2

Mobile unit

for the optimum cleaning of instruments on 2 levels.



- Equipped with an integrated spray arm
- Magnetic strip for automatic mobile unit recognition
- The vertical clearance for both levels is different

EAN: 4002514816449 / Material number: 07765790 / Old Material Number: 69597503D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	6
SPRI mesh trays [number]	4
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in in. (mm)	16 13/16 (427)
External dimensions, net width in in. (mm)	25 1/4 (640)
External dimensions, net depth in in. (mm)	31 1/8 (790)
External dimensions, gross height in in. (mm)	25 5/8 (650)
External dimensions, gross width in in. (mm)	21 11/16 (550)
External dimensions, gross depth in in. (mm)	25 1/4 (640)
Net weight in lb (kg)	24 (10.80)
Gross weight in lb (kg)	29 (12.80)
Accessory pack	
Magnetic strip for automatic program recognition (incl. magnets) [number]	1