

E 198

Insert

for optimum loading of up to 6 mesh and kidney trays.



- Suitable for 6 mesh or kidney trays
- 7 brackets (6 compartments), height 6 5/16 (160) in (mm) mm, distance 2 11/16 (67) in (mm) mm

EAN: 4002514953373 / Material number: 07984040 / Old Material Number: 68519801D

Product category	
Thermal disinfectors	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for mesh trays	•
Insert for kidney dishes / basins	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	6
Distance between rails in mm	2 11/16 (67)
Rail height in mm	6 5/16 (160)
Application area	
Reprocessing of dental instruments	•
Capacity	
Diverse mesh trays [number]	6
Kidney dishes [number]	6
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in in. (mm)	6 5/16 (160)
External dimensions, net width in in. (mm)	7 1/8 (180)
External dimensions, net depth in in. (mm)	19 1/2 (495)
External dimensions, gross height in in. (mm)	7 5/16 (185)
External dimensions, gross width in in. (mm)	8 5/16 (210)
External dimensions, gross depth in in. (mm)	20 7/8 (530)
Net weight in lb (kg)	2 (0.65)
Gross weight in lb (kg)	2 (0.90)