


E 197

Insert

For optimum accommodation of different instruments.



- Mesh width on the base 3 mm, sides 1.7 mm, lid 8 mm
- Can be steam sterilized at 121 °C/134 °C

EAN: 4002514953380 / Material number: 07984090 / Old Material Number: 68519701D

Product classification	
Thermal disinfectors	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of dental instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	42
External dimensions, net width in mm	150
External dimensions, net depth in mm	225
External dimensions, gross height in mm	3 (75)
External dimensions, gross width in mm	6 3/4 (170)
External dimensions, gross depth in in (mm)	9 1/2 (240)
Net weight in kg	2 (0.55)
Gross weight in kg	2 (0.62)