


E 807

Insert

for optimum loading of 3 mesh or kidney trays.



- 4 holders (3 compartments), H 165 mm, spacing 2 11/16 (68) in (mm) mm

EAN: 4002514951638 / Material number: 07984220 / Old Material Number: 69580701D

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
Product category	
Insert for mesh trays	•
Insert for kidney dishes / basins	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	3
Distance between rails in mm	2 11/16 (68)
Application area	
Reprocessing of dental instruments	•
Capacity	
Diverse mesh trays [number]	3
Kidney dishes [number]	3
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	144
External dimensions, net width in mm	140
External dimensions, net depth in mm	290
External dimensions, gross height in mm	18 15/16 (480)
External dimensions, gross width in mm	12 5/8 (320)
External dimensions, gross depth in in (mm)	5 1/8 (130)
Net weight in kg	1 (0.31)
Gross weight in kg	2 (0.47)