


E 806/1

Insert

for the optimum loading of 8 half-shell trays.



- 9 holders (8 compartments), W 295, spacing 1 1/2 (37) in (mm) mm
- Max. tray size 290 x 30 mm

EAN: 4002514969718 / Material number: 09013220 / Old Material Number: 69580603D

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert for trays	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	8
Distance between rails in mm	1 1/2 (37)
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	114
External dimensions, net width in mm	305
External dimensions, net depth in mm	348
External dimensions, gross height in mm	5 (125)
External dimensions, gross width in mm	13 (330)
External dimensions, gross depth in in (mm)	19 (480)
Net weight in kg	2 (0.71)
Gross weight in kg	2 (0.84)