

E 417/1

Insert

For the optimum loading of up to 30 ear and nose specula.

EAN: 4002514984605 / Material number: 09033260 / Old Material Number: 69541702D



- suitable for approx. 30 ear and nose specula
- 280 compartments approx. 10 x 10
- The mesh width on the base is 4 mm
- For use in upper and lower basket

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert for ENT instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	280
<b>Application area</b>	
Reprocessing of ENT instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	3 3/16 (80)
External dimensions, net width in in. (mm)	8 11/16 (220)
External dimensions, net depth in in. (mm)	17 3/16 (435)
External dimensions, gross height in in. (mm)	3 3/8 (85)
External dimensions, gross width in in. (mm)	8 7/8 (225)
External dimensions, gross depth in in. (mm)	18 5/16 (465)
Net weight in lb (kg)	3 (1.30)
Gross weight in lb (kg)	4 (1.72)