

ProCare Universal 61 - 6 kg
 Water softener salt, 3 x 2 kg
 For optimum regeneration of the integrated softener.



- Prevents calcium deposits and protects wash items

EAN: 4002515075319 / Material number: 09195780 /
 Old Material Number: 66916102EU1

Product category	
Fresh-water dishwashers	•
ProfiLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Glassware Washers, laboratory	•
Thermal disinfectors	•
Glassware washers, laboratory	•
Washer-disinfectors	•
Applications	
Reactivation of water softeners	•
Dispensing	
pH value at 20°C (68°F)	10,9
Product characteristics	
Medium type	Reactivation salt
Ingredients	Coarse special salt
Prevents limescale and protects load items	•
Dimensions and weight	
External dimensions, gross depth in mm	225.00
Product characteristics	
Medium type	Reactivation salt
Ingredients	Coarse special salt
Prevents limescale and protects load items	•
Dimensions and weight	
External dimensions, net height in in. (mm)	12 5/8 (320)
External dimensions, net width in in. (mm)	4 1/8 (104)
External dimensions, net depth in in. (mm)	8 15/16 (226)
External dimensions, gross height in in. (mm)	12 5/8 (320)
External dimensions, gross width in in. (mm)	3 15/16 (100)
External dimensions, gross depth in in. (mm)	8 7/8 (225)
Load size in g	211.64 (6000)
Gross weight in lb (kg)	0.21 (6)
Compliance, test certificates, and symbols	
GHS / CLP compliant	•