

E 339/1

Insert

for the optimum loading of up to 12 tray dividers and trays.



- Suitable for 12 tray dividers or trays
- Equipped with 13 holders (12 compartments).
- Maximum tray size 290 x 20 mm

EAN: 4002515094860 / Material number: 09233840 / Old Material Number: 69533902D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert for trays	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	12
Capacity	
Trays [quantity]	12
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	115
External dimensions, net width in mm	305
External dimensions, net depth in mm	498
External dimensions, gross height in mm	18 15/16 (480)
External dimensions, gross width in mm	12 5/8 (320)
External dimensions, gross depth in in (mm)	5 1/8 (130)
Net weight in kg	2 (0.80)
Gross weight in kg	3 (1.22)