

APH 332  
 Adapter cable  
 for the connection of W&H Lisa 317, 322, 517, 522.

EAN: 4002515350959 / Material number: 09553890 / Old Material Number: 68812503D

<b>Product category</b>	
Small sterilizers	•
<b>Product category</b>	
Process documentation / operating documentation / communication	•
<b>Application area</b>	
Reprocessing of pharmacy / infusion bottles	•
Reprocessing anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of mesh trays	•
Reprocessing of containers, e.g., soda lime containers	•
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	0 (0) in (mm)
External dimensions, net width in in. (mm)	0 (0) in (mm)
External dimensions, net depth in in. (mm)	0 (0) in (mm)

- In addition to connection cable type 2
- Length 3 m
- For connecting the machine to Segosoft Miele Edition