

APH 333
 Adapter cable
 for the connection of W&H Lisa MB17, MB22, Madrimed.

EAN: 4002515350966 / Material number: 09554370 / Old Material Number: 68812504D

Product category	
Small sterilizers	•
Product category	
Process documentation / operating documentation / communication	•
Product characteristics	
Colour	Grey
Application area	
Reprocessing of pharmacy / infusion bottles	•
Reprocessing anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of mesh trays	•
Reprocessing of containers, e.g., soda lime containers	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in in. (mm)	0 (0) in (mm)
External dimensions, net width in in. (mm)	0 (0) in (mm)
External dimensions, net depth in in. (mm)	0 (0) in (mm)

- In addition to connection cable type 1
- Length 1 m
- For connecting the machine to Segosoft Miele Edition