


A 100

Upper basket

For the optimum loading of injector modules or inserts.



- Upper basket with two module connections
- for holding up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Loading height with module with plastic supports 190 mm
- Clearances: width 489, length 440 mm

EAN: 4002515460610 / Material number: 09862320 / Old Material Number: 69610000D

<b>Product category</b>	
Glassware Washers, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
<b>Product category</b>	
Basket with module connection	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of pipettes / phials	•
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	5 9/16 (141)
External dimensions, net width in in. (mm)	20 11/16 (525)
External dimensions, net depth in in. (mm)	20 13/16 (528)
External dimensions, gross height in in. (mm)	6 3/4 (170)
External dimensions, gross width in in. (mm)	21 7/8 (555)
External dimensions, gross depth in in (mm)	22 1/4 (565)
Net weight in lb (kg)	7 (2.80)
Gross weight in lb (kg)	8 (3.49)