


A 101

Upper basket

Height adjustable, for the optimum loading of inserts.



- Vertical clearance 160 +/- 30 mm, width 475 mm, depth 443 mm
- With integrated spray arm

EAN: 4002515460627 / Material number: 09862340 / Old Material Number: 69610100D

Product category	
Fresh-water dishwashers	•
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8063	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
Product category	
Upper basket	•
Carrier / basket carrier	•
Upper basket	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Reprocessing of different utensils	•
Reprocessing of air regulators	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	8 1/8 (206)
External dimensions, net width in in. (mm)	20 13/16 (528)
External dimensions, net depth in in. (mm)	20 3/4 (527)
External dimensions, gross height in in. (mm)	9 1/16 (230)
External dimensions, gross width in in. (mm)	21 7/8 (555)
External dimensions, gross depth in in. (mm)	22 7/8 (580)
Net weight in lb (kg)	7 (3.00)
Gross weight in lb (kg)	9 (3.79)