


A 103

Upper basket

For the optimum loading of inserts.



- Vertical clearance: 3 3/4 (95) in (mm), width 505 mm, depth 472 mm
- Particularly suitable for the reprocessing of DIN mesh trays
- Also suitable for further applications with low height items
- Can be used in combination with the A 202 mobile unit
- With integrated spray arm

EAN: 4002515460641 / Material number: 09862370 / Old Material Number: 69610300D

Product category	
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
Product category	
Carrier / basket carrier	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Maximum vertical clearance in mm	3 3/4 (95)
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of OR shoes	•
Reprocessing of surgical instruments	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
Capacity	
DIN mesh trays [number]	2
SPRI mesh trays [number]	1
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	5 1/4 (133)
External dimensions, net width in in. (mm)	20 13/16 (528)
External dimensions, net depth in in. (mm)	20 13/16 (528)
External dimensions, gross height in in. (mm)	6 3/4 (170)
External dimensions, gross width in in. (mm)	21 1/2 (545)
External dimensions, gross depth in in (mm)	22 1/16 (560)
Net weight in lb (kg)	8 (3.50)
Gross weight in lb (kg)	11 (4.57)