


A 104

Upper basket

For the optimum loading of instruments for MIS, arthroscopy and urology.



- Right side for holding 11 hollow instruments
- Additional Luer lock connectors for hollow instruments
- With spray arm: Mesh trays can be placed on the lower carrier

EAN: 4002515460658 / Material number: 09862380 / Old Material Number: 69610400D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Basket for MIS instruments	•
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Capacity	
DIN mesh trays [number]	1
Lumened instruments [number]	11
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	13 1/8 (332)
External dimensions, net width in in. (mm)	20 13/16 (528)
External dimensions, net depth in in. (mm)	20 3/4 (526)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	9 (4.00)
Gross weight in lb (kg)	15 (6.50)