


A 200
Mobile unit
for the optimum loading of 38 pipettes.



- for 10 100 ml pipettes in row 1, distance 20 mm
- for 14 25 ml pipettes in a 2nd row, distance 26 mm
- for 14 10 ml pipettes in a 3. row, distance 26 mm

EAN: 4002515460696 / Material number: 09862420 / Old Material Number: 69620000D

Product category	
Glassware Washers, laboratory	•
Machine affiliation	
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of pipettes / phials	•
Capacity	
Pipettes [number]	38
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	15 11/16 (397)
External dimensions, net width in in. (mm)	20 7/8 (529)
External dimensions, net depth in in. (mm)	21 1/2 (546)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	10 (4.40)
Gross weight in lb (kg)	15 (6.40)