


A 201

Mobile unit

For the optimum loading of anaesthetic sets.



- 8 nozzles with spring holders for breathing tubes up to 1.5 m
- Storage on spiral
- 15 injector nozzles for intubation material

EAN: 4002515460702 / Material number: 09862430 / Old Material Number: 69620100D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Mobile unit for anaesthesia instruments	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Capacity</b>	
AN sets [number]	3
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	17 1/8 (434)
External dimensions, net width in in. (mm)	20 15/16 (531)
External dimensions, net depth in in. (mm)	21 1/2 (546)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in lb (kg)	20 (9.00)
Gross weight in lb (kg)	24 (10.50)
<b>Accessory pack</b>	
E 430/1 1/3 mesh tray [number]	1
E 431 injector nozzle for concertina hose, Ø 8 x 193 mm [number]	1
E 432 holder for concertina breathing tubes [number]	1
E 433 holder for silicone breathing tubes [number]	2
E 434 holder for paediatric breathing tubes [number]	1
E 466 injector nozzle 8 x 333 mm, for breathing bags [number]	6
E 496 injector nozzle 4 x 120 mm [number]	15