


A 204
Mobile unit
for the optimum loading of Micro instruments.



- Level 1: suitable for inserts
- Level 2 right: 22 connections for hollow instruments
- Level 2 left: 11 Luer lock adapters, male, horizontal Hose
- With integrated spray arm
- Tubular filter A 800 can be used in both injector rails

EAN: 4002515460733 / Material number: 09862460 / Old Material Number: 69620400D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Mobile unit for ophthalmology instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Option to use a central filter	•
Application area	
Reprocessing of ophthalmology instruments	•
Capacity	
Surgery sets for micro instruments [number]	2-4
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	14 3/8 (365)
External dimensions, net width in in. (mm)	20 7/8 (529)
External dimensions, net depth in in. (mm)	21 1/8 (535)
External dimensions, gross height in in. (mm)	20 7/8 (530)
External dimensions, gross width in in. (mm)	22 1/16 (560)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in lb (kg)	26 (11.50)
Gross weight in lb (kg)	28 (12.50)
Accessory pack	
E 476 holder [number]	1
E 477 stopper [number]	1