


A 304

Module

For the optimum loading of up to 98 vials.



- Suitable for 98 centrifuge tubes, phials, test tubes, for example
- The loading heights differ depending on the positioning
- 1 unit can be inserted into each of the baskets A 100, A 150, and APLW 106

EAN: 4002515460832 / Material number: 09862570 / Old Material Number: 69630400D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert / module for test tubes	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of pipettes / phials	•
Reprocessing of test tubes	•
Capacity	
Vials [number]	98
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	5 1/8 (130)
External dimensions, net width in in. (mm)	8 3/4 (222)
External dimensions, net depth in in. (mm)	18 9/16 (471)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	6 (2.30)
Gross weight in lb (kg)	8 (3.30)