

A 801  
 Adapter  
 for use of AUF 1 and AUF 2 in connection with A 105/1.



EAN: 4002515460931 / Material number: 09862940 / Old Material Number: 69680100D

<b>Product classification</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
<b>Product category</b>	
Holder and / or adapter for handpieces	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	50
External dimensions, net width in mm	50
External dimensions, net depth in mm	50
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	3 3/16 (80)
External dimensions, gross depth in in (mm)	4 3/4 (120)
Net weight in kg	1 (0.01)
Gross weight in kg	1 (0.01)