

A 802

Jet

for the optimum cleaning of the dispenser compartment.



- For use in laboratory glassware washer with powder dispenser
- For use with injector modules
- Suitable for removing powder residues from dispenser compartment

EAN: 4002515460979 / Material number: 09863070 / Old Material Number: 69680200D

Product classification	
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8504	•
PG 8583	•
Product category	
Other accessories	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	30
External dimensions, net width in mm	187
External dimensions, net depth in mm	15
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	3 3/16 (80)
External dimensions, gross depth in in (mm)	4 3/4 (120)
Net weight in kg	1 (0.02)
Gross weight in kg	1 (0.02)