

A 803  
Holder  
for the optimum loading of transmission instruments.

EAN: 4002515460993 / Material number: 09863100 / Old Material Number: 69680300D



- For use in A 105/1 upper basket along with A 800 tubular filter
- Suitable for numerous contra-angles & hand pieces as well as turbines

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
<b>Product category</b>	
Holder	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	1 5/16 (32) in (mm)
External dimensions, net width in mm	1 5/16 (32) in (mm)
External dimensions, net depth in mm	1 5/8 (40) in (mm)
External dimensions, gross height in mm	3 3/16 (80)
External dimensions, gross width in mm	3 3/4 (95)
External dimensions, gross depth in in (mm)	4 3/8 (110)
Net weight in kg	1 (0.05)
Gross weight in kg	1 (0.09)