

A 800

Tubular filter

For the safe and economic reprocessing of hollow instruments.



- Filter in ply mesh stainless steel mesh for wash carts
- Mesh size 63 µm
- Reprocessable with cleaning brush A 804

EAN: 4002515538081 / Material number: 09960680 / Old Material Number: 68680000D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Filter	•
Product characteristics	
Material	Stainless steel
Dimensions and weight	
External dimensions, net height in mm	0 15/16 (23) in (mm)
External dimensions, net width in mm	0 15/16 (23) in (mm)
External dimensions, net depth in mm	16 (405) in (mm)
External dimensions, gross height in mm	1 13/16 (45)
External dimensions, gross width in mm	8 11/16 (220)
External dimensions, gross depth in in (mm)	17 9/16 (445)
Diameter in in (mm)	0 15/16 (23)
Net weight in kg	1 (0.36)
Gross weight in kg	2 (0.88)