

A 804

Cleaning brush

For the optimum cleaning of the tubular filter A 800.



- Suitable for cleaning tubular filter A 800
- Twisted brushes with spiral trim
- End of handle with eye
- Length 500 mm, brush length 100 mm

EAN: 4002515527238 / Material number: 09960720 / Old Material Number: 68680400D

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Filter	•
Other accessories	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	500
External dimensions, gross height in mm	1 13/16 (45)
External dimensions, gross width in mm	2 13/16 (70)
External dimensions, gross depth in in (mm)	24 1/16 (610)
Diameter in in (mm)	0 13/16 (20)
Net weight in kg	1 (0.07)
Gross weight in kg	1 (0.16)