

PG 8596

Dispensing cabinet

For process chemicals and dispenser modules

EAN: 4002515558324 / Material number: 10087340 / Old Material Number: 69859602D



- 2 removable drawers for 3 x 5 l canisters each
- Internal dimensions H 285-690 , W 250, D 425-555 mm
- Each level 2 holding options per level for dispenser modules
- Drip tray facilitates better residual emptying of the canister

Product category	
Glassware Washers, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Dispenser system	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in in. (mm)	32 7/8 (835)
External dimensions, net width in in. (mm)	11 13/16 (300)
External dimensions, net depth in in. (mm)	23 5/8 (600)
External dimensions, gross height in in. (mm)	35 7/16 (900)
External dimensions, gross width in in. (mm)	15 3/4 (400)
External dimensions, gross depth in in (mm)	29 3/16 (740)
Net weight in lb (kg)	59 (26.50)
Gross weight in lb (kg)	67 (30.00)