

UG 30-60/60-85

Box base

to raise the machine, 60 cm deep, 60 cm wide.



- Suitable for all machines with a width and depth of 60 cm

EAN: 4002515558331 / Material number: 10087350 / Old Material Number: 69306001D

Product category	
Glassware Washers, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8504	•
PG 8583	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Plinth / base	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	12 1/2 (317)
External dimensions, net width in in. (mm)	23 5/8 (600)
External dimensions, net depth in in. (mm)	23 5/8 (600)
External dimensions, gross height in in. (mm)	12 1/4 (310)
External dimensions, gross width in in. (mm)	24 1/16 (610)
External dimensions, gross depth in in (mm)	24 1/16 (610)
Net weight in lb (kg)	22 (9.70)
Gross weight in lb (kg)	25 (11.10)