

MC/2

Transfer cart

height adjustable, for loading & unloading washer-&disinfectors.



- Equipped with 4 castors, 2 of which can be locked
- 2 loading racks (trough shaped) as well as handle and cover plate

EAN: 4002515588628 / Material number: 10138690 / Old Material Number: 69200109D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Storage / transport system	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	39 1/2 (1003)
External dimensions, net width in in. (mm)	23 5/8 (600)
External dimensions, net depth in in. (mm)	31 3/4 (805)
External dimensions, gross height in in. (mm)	44 1/16 (1119)
External dimensions, gross width in in. (mm)	32 7/8 (835)
External dimensions, gross depth in in. (mm)	50 7/16 (1280)
Net weight in lb (kg)	76 (34.20)
Gross weight in lb (kg)	116 (52.40)