


UG 30-90/70-85

Box base

to raise the machine, 70 cm deep, 90 cm wide.



- Suitable for all machines, width 90, depth 70 cm

EAN: 4002515622506 / Material number: 10210600 / Old Material Number: 69309002D

Product category	
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8582 CD	•
PG 8583 CD	•
Product category	
Plinth / base	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	12 1/2 (317)
External dimensions, net width in in. (mm)	35 7/16 (900)
External dimensions, net depth in in. (mm)	27 9/16 (700)
External dimensions, gross height in in. (mm)	12 1/4 (310)
External dimensions, gross width in in. (mm)	24 1/16 (610)
External dimensions, gross depth in in. (mm)	24 1/16 (610)
Net weight in lb (kg)	35 (15.60)
Gross weight in lb (kg)	38 (17.00)