

ProCare Lab 11 AP - 10 kg  
 Powder cleaning detergent, alkaline, 10 kg  
 For optimum reprocessing of laboratory utensils.



- Free from phosphorous compounds, oxidising effect
- Specially developed for Miele Professional
- Suitable for all types of materials
- PH value 11.5-11.7

EAN: 4002515650868 / Material number: 10266540 / Old Material Number: 62922101EU1

Product classification	
Glassware Washers, laboratory	•
Applications	
Removal of inorganic and organic soiling	Good
Dispensing recommendation	
pH value of solution	11.5-11.7
Recommended water range in mmol/l	0-1.1
Main applications	
Water laboratories	•
Serological institutes	•
Phosphate industry	•
Product characteristics	
Medium type	Alkaline cleaning agent
Aggregate form	Powder
Ingredients	Oxygen bleach, Tensides
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for porcelain / ceramics	•
Suitable for laboratory glassware	•
Not suitable for anodized aluminum	•
Does not contain phosphates	•
Not suitable for non-anodized aluminum	•
Not suitable for light and non-ferrous metals	•
Dimensions and weight	
External dimensions, net height in mm	10 9/16 (268)
External dimensions, net width in mm	10 7/16 (265)
External dimensions, net depth in mm	10 7/16 (265)
External dimensions, gross height in mm	11 1/16 (280)
External dimensions, gross width in mm	10 11/16 (270)
External dimensions, gross depth in mm	10 11/16 (270)
Load size in g	352.73 (10000)
Gross weight in kg	0.35 (10)
Compliance, certs / symbols	
GHS / CLP compliant	•