

ProCare Lab 30 C - 5 l

Neutralizing agent, acidic, 5 l

For optimum neutralization with citric acid based neutralizing agent.



- Free from tensides, phosphorus compounds, and oxidation agents
- Specially developed for Miele Professional
- Suitable for all types of materials
- pH value 2.6-3.0
- To be used as neutralising agent or acid pre-cleaner

EAN: 4002515660744 / Material number: 10266670 / Old Material Number: 62927101EU1

Product classification	
Glassware Washers, laboratory	•
Applications	
Removal of tarnishing, limescale, urine scale, and other acid-soluble impurities	Good
Dispensing recommendation	
Dispensing recommendation per rinse	1-4 ml/l (0.1-0.4%)
pH value of solution	2.6-3.0
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Industrial laboratories	•
Phosphate industry	•
Laboratories for the pharmaceutical industry	•
Product characteristics	
Medium type	Acidic neutralising agent
Aggregate form	Liquid
Ingredients	Organic acids
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramics	•
Not suitable for anodized aluminum	•
Does not contain phosphates	•
Does not contain surfactants	•
Very gentle on materials	•
Acidic pre-cleaning	•
Suitable for sensitive materials	•
Neutralization of alkaline carryover	•
Not suitable for chrome-plated parts	•
Not suitable for nickel-plated parts	•
Not suitable for acidic pre-cleaning of non-anodized aluminum	•
Not suitable for light and non-ferrous metals	•
Suitable for neutralization of non-anodized aluminum	•
Dimensions and weight	
External dimensions, net height in mm	8 9/16 (217)
External dimensions, net width in mm	5 5/8 (142)
External dimensions, net depth in mm	9 13/16 (248)
External dimensions, gross height in mm	9 7/8 (250)
External dimensions, gross width in mm	7 1/2 (190)
External dimensions, gross depth in mm	5 15/16 (150)
Gross weight in kg	0.21 (6)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)

ProCare Lab 30 C - 5 l

Neutralizing agent, acidic, 5 l

For optimum neutralization with citric acid based neutralizing agent.

EAN: 4002515660744 / Material number: 10266670 / Old Material Number: 62927101EU1

Compliance, certs / symbols

GHS / CLP compliant

