

ProCare Lab 11 MA - 10 kg

Powder detergent, mildly alkaline, 10 kg

For optimum reprocessing of laboratory utensils.



- Foam damping effect, free from oxidants
- Suitable for the reprocessing phases in washer-disinfectors
- Suitable for all types of materials
- PH value 11.7-11.9

EAN: 4002515650851 / Material number: 10266720 / Old Material Number: 62921101EU1

Product classification	
Glassware Washers, laboratory	•
Applications	
Removal of blood	Good
Removal of protein	Good
Removal of organic residues	Good
Removal of food residues	Good
Dispensing recommendation	
pH value of solution	11.7-11.9
Main applications	
Medical laboratories	•
Biological laboratories	•
Chemical laboratories	•
Product characteristics	
Medium type	Mildly alkaline Cleaning agent
Aggregate form	Powder
Ingredients	
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramics	•
Limited suitability for non-ferrous metals	•
Not suitable for anodized aluminum	•
Reduces foam	•
Does not contain oxidants	•
Limited suitability for non-anodized aluminum	•
Dimensions and weight	
External dimensions, net height in mm	10 9/16 (268)
External dimensions, net width in mm	10 7/16 (265)
External dimensions, net depth in mm	10 7/16 (265)
External dimensions, gross height in mm	11 1/16 (280)
External dimensions, gross width in mm	10 11/16 (270)
External dimensions, gross depth in mm	10 11/16 (270)
Load size in g	352.73 (10000)
Gross weight in kg	0.35 (10)
Compliance, certs / symbols	
GHS / CLP compliant	•