

A 810

Cover net

For optimum securing of delicate items.



- Stainl. steel frames w/ plastic netting, dimensions 445 x 215 mm
- Particularly resistant and durable

EAN: 4002515676783 / Material number: 10283650 / Old Material Number: 69501700D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Cover	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, White
Application area	
Reprocessing of surgical instruments	•
Dimensions and weight	
External dimensions, net height in in. (mm)	0 (0)
External dimensions, net width in in. (mm)	0 (0)
External dimensions, net depth in in. (mm)	0 (0)
External dimensions, gross height in in. (mm)	0 13/16 (20)
External dimensions, gross width in in. (mm)	9 5/16 (235)
External dimensions, gross depth in in (mm)	18 9/16 (470)
Net weight in lb (kg)	1 (0.44)
Gross weight in lb (kg)	2 (0.57)