

A 810/1

Cover net

For optimum securing of delicate items.



- Stainl. steel frames w/ plastic netting, dimensions 206 x 215 mm
- Particularly resistant and durable

EAN: 4002515671764 / Material number: 10283670 / Old Material Number: 69501800D

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Cover	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, White
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	0 7/8 (20)
External dimensions, gross width in mm	9 3/8 (235)
External dimensions, gross depth in in (mm)	9 1/4 (232)
Net weight in kg	1 (0.28)
Gross weight in kg	1 (0.35)