


A 811/1  
Replacement net  
For replacing the net on an A 810/1



- Easy replacement of the netting
- Particularly resistant and durable

EAN: 4002515676493 / Material number: 10283690 / Old Material Number: 69502000D

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Cover	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, White
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	0 7/8 (20)
External dimensions, gross width in mm	9 1/8 (230)
External dimensions, gross depth in in (mm)	9 1/8 (230)
Net weight in kg	1 (0.02)
Gross weight in kg	1 (0.07)