

ProCare Shine 40 - 5 l

Rinse aid, 5 l

For rinsing dishes, cutlery and glasses.

EAN: 4002515692691 / Material number: 10328030 / Old Material Number: 61915201EU1



- Weakly acidic, contains citric acid
- For optimum rinsing results
- For use in commercial dishwashers
- Liquid concentrate, recommendation: 0.02-0.05 % (0.2-0.5 ml/l)
- Specially developed for Miele Professional

Product classification	
Fresh-water dishwashers	•
ProfiLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Dispensing recommendation	
Dispensing recommendation per rinse	0.2-0.5 ml/l (0.02-0.05%)
Large opening for suction lance	•
pH value of solution	6.5-7.0
Recommended water range in mmol/l	0-0.53
Product characteristics	
Medium type	Rinse aid
Aggregate form	Liquid concentrate
Ingredients	Sequestrant (w/o phosphate), Tensides, Citric acid
Suitable for glasses	•
Suitable for cutlery	•
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for porcelain / ceramics	•
Suitable for aluminum	•
Suitable for silver	•
Gentle on decor	•
Prevents limescale and protects load items	•
Improves and shortens the drying process	•
Waste water neutral	•
Dimensions and weight	
External dimensions, net height in mm	9 7/8 (250)
External dimensions, net width in mm	5 15/16 (150)
External dimensions, net depth in mm	7 1/2 (190)
External dimensions, gross height in mm	9 7/8 (250)
External dimensions, gross width in mm	5 15/16 (150)
External dimensions, gross depth in mm	7 1/2 (190)
Gross weight in kg	0.18 (5)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
GHS / CLP compliant	•