

E 454/1

Injector nozzle

for the optimum cleaning of diverse loads, particularly for proctoscopes.



- Suitable for proctoscopes. Trocar sleeves and other wash loads

EAN: 4002515722954 / Material number: 10342490 / Old Material Number: 69745402D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7824	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Injector nozzle	•
Other accessories	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of lumened instruments	•
Capacity	
Lumened instruments [number]	1
Dimensions and weight	
External dimensions, net height in mm	0 3/8 (8) in (mm)
External dimensions, net width in mm	0 3/8 (8) in (mm)
External dimensions, net depth in mm	11 13/16 (300) in (mm)
External dimensions, gross height in mm	2 3/8 (60)
External dimensions, gross width in mm	6 5/16 (160)
External dimensions, gross depth in in (mm)	13 13/16 (350)
Diameter in in (mm)	0 3/8 (8)
Net weight in kg	1 (0.06)
Gross weight in kg	1 (0.07)