

ProCare Shine 11 OB - 6 kg
 Powder detergent, mildly alkaline, 6 kg
 For optimum cleaning of dishes, cutlery, and glasses.

EAN: 4002515802212 / Material number: 10582700 / Old Material Number: 61911202EU1



- Contains active oxygen and phosphates
- Excellent cleaning performance, particularly on coffee/tea stains
- For use in commercial dishwashers
- Dispensing recommendation: 20–25 g per wash cycle or 2–5 g/l
- Specially developed for Miele Professional

Product classification	
Fresh-water dishwashers	•
ProfiLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Applications	
Removal of coffee and tea residues	Very good
Removal of starch/protein	Very good
Removal of dried food residues	Very good
Removal of burnt-on grease	Very good
Removal of grease	Very good
Dispensing recommendation	
Dispensing recommendation per rinse	20–25 g / 2–5 g/l
Intensity of soiling	Strong
pH value of solution	11.0-11.5
Recommended water range in mmol/l	0-0.53
Product characteristics	
Medium type	Mildly alkaline Cleaning agent
Aggregate form	Powder
Ingredients	Alkalis,Oxygen bleach,Phosphate
Suitable for glasses	•
Suitable for cutlery	•
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for porcelain / ceramics	•
Suitable for aluminum	•
Suitable for silver	•
Gentle on decor	•
Protects against corrosion	•
Prevents limescale and protects load items	•
Dimensions and weight	
External dimensions, net height in mm	8 7/16 (214)
External dimensions, net width in mm	9 3/16 (233)
External dimensions, net depth in mm	9 3/16 (233)
External dimensions, gross height in mm	8 7/16 (214)
External dimensions, gross width in mm	9 3/16 (233)
External dimensions, gross depth in mm	9 3/16 (233)
Diameter in mm	9 3/16 (233)
Load size in g	211.64 (6000)
Gross weight in kg	0.21 (6)
Compliance, certs / symbols	
GHS / CLP compliant	•