

A 815
Nozzle

Perfect reprocessing of transmission instruments with an external spray channel.

EAN: 4002515830413 / Material number: 10619690 / Old Material Number: 68681500D



- For use in A 812, A 813, or in A 105/1 upper basket
- Consists of 4 connectors with silicone hoses, L approx. 114 mm

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Transparent
Application area	
Reprocessing of dental instruments	•
Capacity	
Handpieces [number]	4
Dimensions and weight	
External dimensions, net height in in. (mm)	0 7/16 (10)
External dimensions, net width in in. (mm)	0 7/16 (10)
External dimensions, net depth in in. (mm)	4 1/2 (114)
External dimensions, gross height in in. (mm)	3 3/16 (80)
External dimensions, gross width in in. (mm)	3 3/4 (95)
External dimensions, gross depth in in. (mm)	4 3/8 (110)
Diameter in in (mm)	0 7/16 (10)
Net weight in lb (kg)	1 (0.05)
Gross weight in lb (kg)	1 (0.12)