

E 473/2

Filter

with lid, suitable for micro components.



- For mounting
- With smooth cover closure

EAN: 4002515811962 / Material number: 10623720 / Old Material Number: 69547303D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
Product category	
Other inserts	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Dimensions and weight	
External dimensions, net height in mm	3 1/2 (88) in (mm)
External dimensions, net width in mm	2 5/8 (66) in (mm)
External dimensions, net depth in mm	2 5/8 (66) in (mm)
External dimensions, gross height in mm	3 3/16 (80)
External dimensions, gross width in mm	3 3/4 (95)
External dimensions, gross depth in in (mm)	4 3/8 (110)
Net weight in kg	1 (0.10)
Gross weight in kg	1 (0.14)