

A 207

Mobile unit

For the optimum loading of Micro instruments.



- Efficient reprocessing: fast loading and high throughput
- Top level: up to 8 trays for hollow instrument loading
- Lower level: 8 horizontal trays for further ophthal. instruments

EAN: 4002515815601 / Material number: 10632400 / Old Material Number: 69620700D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Mobile unit for ophthalmology instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Option to use a central filter	•
Application area	
Reprocessing of ophthalmology instruments	•
Capacity	
Trays [quantity]	16
Surgery sets for micro instruments [number]	8
Dimensions and weight	
External dimensions, net height in in. (mm)	17 1/4 (438)
External dimensions, net width in in. (mm)	20 7/8 (530)
External dimensions, net depth in in. (mm)	21 3/8 (542)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in lb (kg)	10 (4.40)
Gross weight in lb (kg)	15 (6.40)
Accessory pack	
A 311 tray holder [quantity]	2
A 800 tubular filter [quantity]	1
A 804 cleaning brush [quantity]	1