

A 310
Insert
for the optimum loading of 8 OR shoes.



- For use in upper baskets
- For loading 8 OR shoes up to size 41
- Can be used in upper baskets A 101 and A 102

EAN: 4002515816349 / Material number: 10632610 / Old Material Number: 69631000D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for OR shoes	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of OR shoes	•
Capacity	
OR shoes [number]	8
Dimensions and weight	
External dimensions, net height in in. (mm)	3 (75)
External dimensions, net width in in. (mm)	16 (405)
External dimensions, net depth in in. (mm)	17 3/16 (435)
External dimensions, gross height in in. (mm)	3 15/16 (100)
External dimensions, gross width in in. (mm)	16 3/4 (425)
External dimensions, gross depth in in (mm)	17 3/4 (450)
Net weight in lb (kg)	2 (0.50)
Gross weight in lb (kg)	2 (0.83)