

A 832

Injector nozzle

for the optimum cleaning of a variety of load items.

EAN: 4002515817711 / Material number: 10635440 / Old Material Number: 69745202D



- Suitable for injector carts, with screw connection
- Optimum cleaning results thanks to lateral drill holes

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Injector nozzle	•
Application area	
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net depth in in. (mm)	2 3/8 (60)
External dimensions, gross height in in. (mm)	0 1/2 (12)
External dimensions, gross width in in. (mm)	3 3/16 (80)
External dimensions, gross depth in in. (mm)	4 3/4 (120)
Diameter in in (mm)	0 1/8 (3)
Net weight in lb (kg)	1 (0.00)
Gross weight in lb (kg)	1 (0.01)