

A 834

Injector nozzle

for the optimum cleaning of diverse load items, particularly for proctoscopes.



- Suitable for proctoscopes, trocar sleeves, and other load items
- Equipped with a spring for height adjustment
- For optimum cleaning results thanks to lateral drill holes

EAN: 4002515940709 / Material number: 10635650 /
Old Material Number: 69745403D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Injector nozzle	•
Application area	
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net width in in. (mm)	0 3/8 (8)
External dimensions, net depth in in. (mm)	6 (152)
External dimensions, gross height in in. (mm)	0 1/2 (12)
External dimensions, gross width in in. (mm)	4 3/4 (120)
External dimensions, gross depth in in. (mm)	6 3/4 (170)
Net weight in lb (kg)	1 (0.03)
Gross weight in lb (kg)	1 (0.03)