

A 817

Irrigation sleeves

for optimum cleaning of ENT instruments, for example.

EAN: 4002515880593 / Material number: 10706260 / Old Material Number: 69681700D



- Suitable for injector carts and baskets, screw connection
- Incl. cap and mounting clip

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8534	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	4 3/4 (120)
External dimensions, gross depth in in (mm)	6 3/4 (170)
Diameter in in (mm)	1 (24)
Net weight in kg	1 (0.02)
Gross weight in kg	1 (0.02)