

A 818

Irrigation sleeve cover

for the optimum securing of instruments in irrigation sleeves.



- Packaging unit with 12 caps
- Opening can be adapted to different applications

EAN: 4002516100119 / Material number: 10706620 / Old Material Number: 69681800D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Injector nozzle	•
Application area	
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Product characteristics	
Material	Plastic
Colour	White
Dimensions and weight	
External dimensions, net height in in. (mm)	0 13/16 (20)
External dimensions, net width in in. (mm)	0 13/16 (20)
External dimensions, net depth in in. (mm)	0 13/16 (20)
External dimensions, gross height in in. (mm)	0 1/2 (12)
External dimensions, gross width in in. (mm)	3 3/16 (80)
External dimensions, gross depth in in. (mm)	4 3/4 (120)
Diameter in in (mm)	0 13/16 (20)
Net weight in lb (kg)	0 (0.04)
Gross weight in lb (kg)	0 (0.04)