

A 819

Irrigation sleeves

for optimum cleaning of lumened instruments.

EAN: 4002515990278 / Material number: 10721750 / Old Material Number: 68681900D



- E.g., for magnetostrictive ultrasound scaler tips, e.g., Cavitron
- Int. dimension 108 mm
- Including cap with an opening of Ø 10 mm
- With mounting clip

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8528	•
G 7831	•
PG 8536	•
PG 8581	•
PG 8591	•
Product category	
Holder	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of dental instruments	•
Capacity	
Ultrasonic tips [number]	1
Dimensions and weight	
External dimensions, net height in in. (mm)	0 7/16 (10)
External dimensions, net width in in. (mm)	3 15/16 (100)
External dimensions, net depth in in. (mm)	4 1/2 (114)
External dimensions, gross height in in. (mm)	0 1/2 (12)
External dimensions, gross width in in. (mm)	4 3/4 (120)
External dimensions, gross depth in in. (mm)	6 3/4 (170)
Net weight in lb (kg)	1 (0.02)
Gross weight in lb (kg)	1 (0.02)