

DOS G 80 Flex Peristaltic pump

With level indicator for the automatic dispensing of liquid cleaning agents.



- For canisters up to 340 mm tall (up to canister opening)
- Dispenses up to 55 Milliliter ml per minute
- Adjustable via the dishwasher controls
- Flexible vent valve cap
- 13 7/8 (352) in (mm) mm suction lance for 5 and 10 l media canisters

EAN: 4002516280064 / Material number: 10730850 / Old Material Number: 67747040USA

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Glassware Washers, laboratory	•
Machine affiliation	
PG 8056	•
PG 8056 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8583	•
PG 8593	•
Product category	
External dispenser system	•
Product characteristics	
Material	Plastic
Colour	Grey
Production rate in ml/60 s	55
Length of suction lance in mm	13 7/8 (352)
Length of dispensing hose to suction lance in inches (mm)	1800
Length of dispensing hose from machine to dispenser module in inches (mm)	2850
Length of power supply cable to machine in inches (mm)	2800
Minimum width of canister opening in inches (mm)	31
Maximum width of canister opening in inches (mm)	60
Standard electrical connection	
Number of phases	1
Neutral conductor	N
AC / DC	AC
Voltage in V	200-240
Frequency in Hz	60
Total rated load in kW	0.01
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	4 1/2 (114) in (mm)
External dimensions, net width in mm	4 9/16 (115) in (mm)
External dimensions, net depth in mm	5 1/2 (139) in (mm)
External dimensions, gross height in mm	6 5/16 (160)
External dimensions, gross width in mm	7 13/16 (198)
External dimensions, gross depth in in (mm)	18 7/8 (478)
Net weight in kg	6 (2.27)
Gross weight in kg	6 (2.51)