


APFD 402

Insert

for the optimum loading of medication cup lids.



- Suitable for loading 8 lids for medication cups
- For use in basket insert APFD 401
- Suitable for lids with a diameter of 28 to 35 mm

EAN: 4002516152170 / Material number: 10869850 / Old Material Number: 67540201D


Product category	
Fresh-water dishwashers	•
Machine affiliation	
PG 8080	•
PG 8080 i	•
PG 8080 U	•
PG 8080 U/i	•
PG 8081 i	•
PG 8082 Sci XXL	•
PG 8083 SCVi XXL	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
Product category	
Insert for medication cup lids	•
Product characteristics	
Material	Plastic
Colour	Grey
Capacity	
Medication cup lids [number]	8

APFD 402

Insert

for the optimum loading of medication cup lids.

EAN: 4002516152170 / Material number: 10869850 / Old Material Number: 67540201D



Dimensions and weight	
External dimensions, net height in mm	0 (0) in (mm)
External dimensions, net width in mm	0 (0) in (mm)
External dimensions, net depth in mm	0 (0) in (mm)
External dimensions, gross height in mm	2 3/16 (55)
External dimensions, gross width in mm	7 1/8 (180)
External dimensions, gross depth in in (mm)	9 7/8 (250)
Net weight in kg	1 (0.01)
Gross weight in kg	1 (0.02)